


## Compact laboratory balance KERN PCB



### The standard in the laboratory



**GLP/ISO record keeping** of weighing data, balance adjustment etc. with date, time and identification no. in combination with an appropriate KERN printer. Ideal for monitoring and documenting your processes in accordance with your quality management system



**Plus/minus weighing procedures**  
For example for checking piece weight, production control, for manufacturing packages of equal weight, etc. The deviation of the test items from a defined target weight will be displayed using the respective "+" and "-" symbols.



With the **recipe function** you can weigh the different ingredients of a mixture. As a check, you can also call up the total weight of all the ingredients

# Compact laboratory balance KERN PCB



## Features

- **PRE-TARE function** for manual subtraction of a known container weight, useful for checking fill-levels
- **Freely programmable weighing unit**, e.g. display direct in special units such as length of thread g/m, paper weight g/m<sup>2</sup>, or similar
- **Percentage determination**: makes it possible to store a given weight value (100 %) and to determine deviations from this target value
- **Ring-shaped draught shield** standard, only for models with weighing plate size **A**, weighing space ØxH 90x40 mm

## Technical data

- Backlit LCD display, digit height 15 mm
- Weighing plate dimensions, stainless steel,
  - A** Ø 81 mm
  - B** Ø 105 mm
  - C** WxD 130x130 mm
  - D** WxD 150x170 mm, see enlarged picture
- Weighing plate material
  - A** Plastic, with conductive lacquer
  - B, C, D** Stainless steel
- Optional battery operation, battery 9 V block not standard. Operating time up to 12 h. AUTO-OFF function to preserve the battery, can be switched off
- Overall dimensions (without draught shield) WxDxH 163x245x79 mm
- Permissible ambient temperature 5 °C / 35 °C

## Accessories

- **Protective working cover** over keyboard and housing, standard. Can be re-ordered, scope of delivery: 5 items for models with weighing plate size
  - A** KERN PCB-A02S05
  - B** KERN PCB-A03S05
  - C** KERN PCB-A04S05
  - D** KERN PCB-A05S05
- **Hook for underfloor weighing** of hanging loads, standard, can be reordered, KERN 440-A01
- **Rechargeable battery pack internal**, can be retrofitted, operating time up to 48 h without backlight, charging time approx. 8 h. AUTO-OFF function to preserve the battery, can be switched off, KERN PCB-A01
- **Software Balance Connection**, for details see page 179, KERN SCD-4.0
- **Individual header data**: the free software KERN SHM-01 can be used to set up 4 header lines on the printout for printers YKN-01, 911-013 and YKB-01N
- **RS-232/Ethernet adapter** for connection to an IP-based Ethernet network, for details see page 180, KERN YKI-01
- **Suitable test weights**, also with calibration certificate see page 188
- **Suitable printers** see page 177 ff.

### STANDARD




### OPTION



Model	Weighing range [Max] g	Readout [d] g	Reproducibility g	Linearity g	Min. piece weight [Counting] g/piece	Net weight approx. kg	Weighing plate	Option	
								DAkkS Calibr. Certificate	
<b>KERN</b>								<b>DKD</b>	
<b>PCB 100-3</b>	100	0,001	0,001	± 0,003	0,002	1,1	<b>A</b>	963-127	
<b>PCB 250-3</b>	250	0,001	0,001	± 0,005	0,002	1,1	<b>A</b>	963-127	
<b>PCB 350-3</b>	350	0,001	0,002	± 0,005	0,002	1,1	<b>A</b>	963-127	
<b>PCB 200-2</b>	200	0,01	0,01	± 0,02	0,02	1,1	<b>B</b>	963-127	
<b>PCB 1000-2</b>	1000	0,01	0,01	± 0,03	0,02	1,4	<b>C</b>	963-127	
<b>PCB 2500-2</b>	2500	0,01	0,01	± 0,5	0,02	1,4	<b>C</b>	963-127	
<b>PCB 3500-2</b>	3500	0,01	0,02	± 0,05	0,02	1,4	<b>C</b>	963-127	
<b>PCB 1000-1</b>	1000	0,1	0,1	± 0,2	0,2	1,4	<b>C</b>	963-127	
<b>PCB 2000-1</b>	2000	0,1	0,1	± 0,2	0,2	1,4	<b>C</b>	963-127	
<b>PCB 6000-1</b>	6000	0,1	0,1	± 0,3	0,2	2	<b>D</b>	963-128	
<b>PCB 10000-1</b>	10000	0,1	0,1	± 0,3	0,2	2	<b>D</b>	963-128	
<b>PCB 6000-0</b>	6000	1	1	± 2	2	2	<b>D</b>	963-128	

# KERN Pictograms:

 <b>Internal adjusting:</b> Quick setting up of the balance's accuracy with internal adjusting weight (motordriven).	 <b>Piece counting:</b> Reference quantities selectable. Display can be switched from piece to weight.	 <b>Suspended weighing:</b> Load support with hook on the underside of the balance.
 <b>Adjusting program CAL:</b> For quick setting up of the balance's accuracy. External adjusting weight required.	 <b>Recipe level A:</b> Separate memory for the weight of the tare container and the recipe ingredients (net total).	 <b>Battery operation:</b> Ready for battery operation. The battery type is specified for each device.
 <b>Memory:</b> Balance memory capacity, e.g. for article data, weighing data, tare weights, PLU etc.	 <b>Recipe level B:</b> Internal memory for complete recipes with name and target value of the recipe ingredients. User guidance through display.	 <b>Rechargeable battery pack:</b> Rechargeable set.
 <b>Alibi memory:</b> Electronic archiving of weighing results, complying with the 2009/23/EC standard.	 <b>Recipe level C:</b> Internal memory for complete recipes with name and target value of the recipe ingredients. User guidance through display, adjustment of recipe when dosages are exceeded, multiplier function, barcode.	 <b>Universal mains adapter:</b> with universal input and optional input socket adapters for A) EU, GB B) EU, GB, CH, USA C) EU, GB, CH, USA, AUS
 <b>Data interface RS-232:</b> To connect the balance to a printer, PC or network.	 <b>Totalising level A:</b> The weights of similar items can be added together and the total can be printed out.	 <b>Mains adapter:</b> 230V/50Hz in standard version for EU. On request GB, USA or AUS version available.
 <b>RS-485 data interface:</b> To connect the balance to a printer, PC or other peripherals. High tolerance against electromagnetic disturbance.	 <b>Totalising level C:</b> Internal memory for complete recipes with name and target value of the recipe ingredients. User guidance through display, adjustment of recipe when dosages are exceeded, multiplier function, barcode recognition.	 <b>Power supply:</b> Integrated in balance. 230V/50Hz standard EU. More standards e.g. GB, USA or AUS on request.
 <b>USB data interface:</b> To connect the balance to a printer, PC or other peripherals.	 <b>Weighing principle: Strain gauge</b> Electrical resistor on an elastic deforming body.	 <b>Weighing principle: Tuning fork</b> A resonating body is electromagnetically excited, causing it to oscillate.
 <b>Bluetooth* data interface:</b> To transfer data from the balance to a printer, PC or other peripherals.	 <b>Percentage determination:</b> Determining the deviation in % from the target value (100 %).	 <b>Weighing principle: Electromagnetic force compensation</b> Coil inside a permanent magnet. For the most accurate weighings.
 <b>WLAN data interface:</b> To transfer data from the balance to a printer, PC or other peripherals.	 <b>Weighing units:</b> Can be switched to e.g. non-metric units at the touch of a key. See balance model. Please refer to KERN's website for more details.	 <b>Weighing principle: Single cell technology</b> Advanced version of the force compensation principle with the highest level of precision.
 <b>Control outputs (optocoupler, digital I/O):</b> To connect relays, signal lamps, valves, etc.	 <b>Weighing with tolerance range:</b> Upper and lower limiting values can be programmed individually for e.g. dosing, sorting and portioning.	 <b>Verification possible:</b> The time required for verification is specified in the pictogram.
 <b>Interface for second balance:</b> For direct connection of a second balance.	 <b>Hold function:</b> (Animal weighing program) When the weighing conditions are unstable, a stable weight is calculated as an average value.	 <b>DAKkS calibration possible (DKD):</b> The time required for DAKkS calibration is shown in days in the pictogram.
 <b>Network interface:</b> For connecting the scale to an Ethernet network. With KERN products you can use a universal RS-232/LAN converter.	 <b>Protection against dust and water splashes IPxx:</b> The type of protection is shown in the pictogram.	 <b>Package shipment:</b> The time required for internal shipping preparations is shown in days in the pictogram.
 <b>Wireless data transfer:</b> between the weighing unit and the evaluation unit using an integrated radio module.	 <b>ATEX explosion protection:</b> Suitable for use in hazardous industrial environments, in which there is explosion danger. The ATEX marking is specified for each device.	 <b>Pallet shipment:</b> The time required for internal shipping preparations is shown in days in the pictogram.
 <b>GLP/ISO log:</b> The balance displays the weight, date and time, regardless of a printer connection.	 <b>Stainless steel:</b> The balance is protected against corrosion.	 <b>Warranty:</b> The warranty period is shown in the pictogram.
 <b>GLP/ISO log:</b> With weight, date and time. Only with KERN printers.		

## KERN – Precision is our business

To ensure the high precision of your balance KERN offers you the the appropriate test weight in the international OIML error limit classes E1-M3 from 1 mg - 2000 kg. In combination with a DAKkS calibration certificate the best pre-requisite for proper balance calibration.

The KERN DAKkS calibration laboratory today is one of the most modern and best-equipped DAKkS calibration laboratories for balances, test weights and force-measurement in Europe.

Thanks to the high level of automation, we can carry out DAKkS calibration of

balances, test weights and force-measuring devices 24 hours a day, 7 days a week.

### Range of services:

- DAKkS calibration of balances with a maximum load of up to 50 t
- DAKkS calibration of weights in the range of 1 mg – 2500 kg
- Database supported management of checking equipment and reminder service
- Calibration of force-measuring devices
- DAKkS calibration certificates in the following languages D, GB, F, I, E, NL, PL

## Your KERN specialist dealer:



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

