





## LARGE CAPACITY FREESTANDING PHARMACY FRIDGE GLASS DOOR 1200L WHITE



- 
**+2°C to +8°C Operating Temperature Range**  
 Conforms to MHRA guidelines.
- 
**NEW Advanced Controller**  
 Records temperatures accurate to 0.1°C.
- 
**Fitted Lock as Standard**  
 Maintains security and ensures the unit can only be accessed by approved personnel.
- 
**Fitted Castors**  
 For ease of movement.
- 
**Wire Shelves**  
 Enables optimum air flow inside the unit.

SKU CODE:      PRODUCT CODE:

444440366      PGR1200UK

DESCRIPTION:

Large Capacity Freestanding Pharmacy Fridge Glass Door 1200L White



## FULL SPECIFICATION – PGR1200UK LARGE CAPACITY FREESTANDING PHARMACY FRIDGE GLASS DOOR 1200L WHITE

<b>Description</b>	Freestanding
<b>SKU code</b>	444440366
<b>GDPA code</b>	PGR1200UK
<b>Gross capacity (litres)</b>	1200
<b>Unit external width (mm)</b>	1330
<b>Unit external depth (mm)</b>	840
<b>Unit external height (mm)</b>	1995
<b>Product weight (kg)</b>	210
<b>Packed width (mm)</b>	1480
<b>Packed depth (mm)</b>	990
<b>Packed height (mm)</b>	2195
<b>Packed weight (kg)</b>	230
<b>Country of origin</b>	UK
<b>Internal dimensions (mm) (WxDxH)</b>	1220 x 653 x 1386
<b>Temperature range</b>	+2°C to +8°C
<b>Colour</b>	White/green
<b>External finish</b>	White
<b>Internal finish</b>	Stainless steel
<b>Type of plastic</b>	SRL800 (HIPS/ABS Mix)
<b>Door finish</b>	Painted stainless steel
<b>Door type</b>	Glass
<b>Door style</b>	Hinged
<b>Door open depth</b>	1145
<b>Lock</b>	Yes
<b>Adjustable feet (number &amp; position)</b>	Fitted castors
<b>Temperature display on panel</b>	Yes
<b>Temperature alarm</b>	Yes
<b>Number of shelves</b>	10
<b>Number of shelf positions</b>	44
<b>Shelf size</b>	645 x 523
<b>Shelf material</b>	Plastic coated wire
<b>Interior light</b>	Yes
<b>Fan</b>	Yes
<b>Spark free</b>	No
<b>Auto or manual defrost</b>	Automatic
<b>Controller type</b>	Digital, external
<b>Climate class</b>	N
<b>Refrigerant type</b>	R134a
<b>CFC free</b>	Yes
<b>HFC free</b>	Yes
<b>HCFC free</b>	Yes
<b>Ammonia free</b>	Yes
<b>UK 3 pin plug</b>	Yes
<b>Other plug types available</b>	Yes
<b>Voltage supply</b>	230
<b>Aux output</b>	Yes
<b>Insulation thickness (mm)</b>	60
<b>Insulation type</b>	Cyclopentane
<b>Anti microbial handle</b>	TBC
<b>Optional drawer system available</b>	Yes



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

