



# LABCOLD™

medical & scientific refrigeration



## RAFG44042 Glass Door Advanced Laboratory Refrigerator

Part Number	RAFG44042
Capacity (Litres)	1245
Style	Freestanding
Exterior Dimensions (HxWxD)	1980 x 1440 x 820mm
Energy Consumption	8.2 kWh/24h
Operating Temperature	+4°C
Temperature Set Point Range	+2°C to +8°C
Interior Construction	Stainless Steel
Exterior Construction	Powder Coated Stainless Steel
Shelves	8
Shelf dimensions (W x D)	530 x 650mm
Door Hinging	Opens from centre outwards
Door Style	Glass
Internal Light	✓
Front Lockable Castors	✓
2 Door Locks	✓
Access Port	Optional Extra
Digital Temperature Control	✓
Digital Temperature Display	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Door Open Alarm	✓
Maximum Temperature Display	✓
Minimum Temperature Display	✓
Auto Defrost	✓
Fan Assisted Cooling	✓
Environmentally Friendly Insulation	✓
Eco Friendly Refrigerant	✓
Ammonia Free	✓

Robustly constructed from stainless steel this refrigerator from the **Labcold** Advanced Range is designed to provide close temperature control for temperature sensitive items.

In addition, remote alarm contacts can be fitted as well as access ports.



The air is continually circulated internally to ensure that where ever product is placed inside the cabinet, it stays within range. There is no internal 'step' and the interior is curved making it easy to clean and stock. The 8 heavy duty shelves are fully adjustable and more shelves can be purchased as an optional extra. The robust doors are independently lockable and for increased security and convenience, are self closing.

Insulated with Envirofoam, which is derived from plant rather than petrochemical oils, the refrigerant is also ammonia, HCFC and CFC free.

The cabinet is fitted with heavy duty, lockable castors making it easy to manoeuvre and clean behind which means this Advanced Laboratory refrigerator is ideal for demanding environments.





**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

