



## RLDG1010 Glass Door Pharmacy & Vaccine Refrigerator

Part Number	RLDG1010
Capacity (Litres)	300
Style	Free standing
Exterior Dimensions (HxWxD)	1500 x 600 x 660mm
Interior Dimensions (HxWxD)	1370 x 515 x 465mm
Energy Consumption	0.67 kWh/24h
Operating Temperature	+5°C
Interior Construction	Painted Corrosion proof Aluminium
Exterior Construction	Powder Coated Steel
Shelves	6
Internal Drawers Available	✓
Door Hinging	Right hand only
Door Style	Glass
Door Lock	✓
Access Port	✓
Digital Temperature Control	✓
Digital Temperature Display	✓
Internal Light	✓
High/Low Temperature alarm	✓
Minimum/Maximum Temperature Recording	✓
Up to 13 months temperature data storage	✓
Door Open alarm	✓
Mains Failure alarm	✓
Alarm Battery Back-Up	✓
Off cycle auto defrost	✓
Fan Assisted Cooling	✓
Remote Alarm Contacts	✓
Eco Friendly Refrigerant	✓
Ammonia Free	✓
DH Green Book, RPS/ GPhC, NaTHNaC and RCVS Guidelines compliant	✓

Powered by the unique IntelliCold™ Controller, this pharmacy refrigerator has superior temperature uniformity and energy consumption, plus extra functions designed to store temperature sensitive pharmaceuticals in optimum conditions and assist busy professionals with regulatory requirements.



The simple to use four button interface is clearly labelled and allows minimum/maximum temperatures to be displayed at the touch of a button. It also shows clearly which alarm is active and why, making remedial action easy. A clear LED display shows the internal temperature of the fridge to one decimal place.

For objective and accurate record keeping, once the warranty is registered, a complimentary SD card will be sent which can be used to download all the temperature data stored by the controller. This can also be used to extract short term data if, for example, there is a suspected temperature incursion; data can be downloaded and used to prove the viability of the contents. Extra reassurance is provided by controller battery back up which can record and alarm in the event of a mains power failure.



The RLDG1010 has a toughened safety glass door and the external light switch on the IntelliCold™ Controller means that the light can operate independent of the door opening.

We are confident that this pharmacy refrigerator will provide you with many years of fault free service. However, for extra peace of mind your RLDG1010 comes complete with a FREE 5 year parts and labour warranty!\*



\* Terms and conditions apply, please visit the web site



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

