



Battery and Filter not included in warranty



**Soft-Coated Triggers**  
with direction indicators

**Ergonomically Balanced**  
in the hand with a loaded pipette

**Universal Grip**  
for right- or left- handed users



**Adapter**  
for small pipettes

### Easy Identification

5 different colored nose cones help differentiate department, task or user; as well as a labeling area on the side of the unit



### Magnetic & Screw Wall-Mount



### 3 Modes

Bright LEDs indicate 3 operation modes: high, low, & gravity



## RF3000®

Packed with features for comfortable and accurate pipetting

The powerful pump offers rapid, yet sensitive aspiration and dispensing, has adjustable valves and is microprocessor controlled for constant reliability and accuracy

- Long-life li-ion battery, with quick charge—up to 20 hours use
- Toggle operation by pressing the large mode button
- An adapter can be fitted to the nose cone to stabilize small volume and thinner pipettes

- New USB connector ensures the appropriate charger is used with the unit
- Low-battery indicator
- Can be recharged when in use



Nose cone assembly only



**Storage Stand**  
doubles as a charging station

### Specifications

<b>Power</b>	Micro-USB connector Rechargeable li-Ion battery (replaceable)
<b>Filter</b>	PTFE/25mm dia/0.45 µm filter
<b>Pipettes</b>	1–100 mL, glass and plastic
<b>Pipette Nozzle</b>	Polypropylene



### Package Includes

- RF3000® Pipette Controller
- Charging Stand and Wall-Mountable/Magnetic Holder
- Nozzle Adapter for use with Small Pipettes
- Li-Ion Rechargeable Battery
- Five Nose Cones: Black, Yellow, Red, Grey and Blue
- 5VDC Micro USB Plug
- 6 ft (183 cm) Cable
- 1 Low Voltage, Double Insulated Power Adapter (100-240 VAC, 50/60 Hz) with 4 Interchangeable Plugs



Item No.	Volume Range mL	L x W x H in	L x W x H cm
HS3000	1–100	9 x 8.5 x 5	22.9 x 21.6 x 12.7

### Replacement parts

HS5000RFFIL	RF3000® Replacement Filters, Pack of 5		
120669	RF3000® Universal Power Adaptor 5VDC Micro USB Plug		
HS300502	RF3000® Li-Ion Rechargeable Battery, One year warranty on battery included		
120487	RF3000® Nose Cones Assorted* ■ ■ ■ ■ ■ Set of 5		
HS5000RFRFA	RF3000® Replacement Filter Adapter		
HS5000RFRP	RF3000® Replacement Pipette Adapter		

\*Assorted colors (Black, Blue, Grey, Red and Yellow)



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

