

# A Precise Electronic Counter and Permanent Sharpie® Marker Combined in One



## eCount™ Colony Counter

Ideal for rapid, accurate counts

All-in-one eCount combines a multi-function electronic counter with a Sharpie® pen for marking Petri dishes to prevent missing or double-counting colonies. Light, ergonomic body is balanced for easy manipulation.

Touch controls for multiple counts, sums, counting up or down, and more! Counter automatically beeps with each count—or select the LED flash mode for a visual alert of each count. Press the lock key to write with Sharpie pen without counting—makes writing count totals simple.

Sharpie pen provides highly visible marks on both plastic and glass dishes. Counter accepts both fine and extra fine point Sharpie pens. Easily change colors for separate counts. Includes one fine-point Sharpie pen and weighted stand.

Part Number - HSI5092A

Manufactured by:  **Heathrow Scientific** LLC

Patent Pending



Included



In use

### eCount Features:

-  Save battery life.
-  Verify each count with either beep or LED flash.
-  First count up, then count down to verify accuracy.
-  Display the total of all counts (up to 32 separate counts).
-  Keep up to 32 separate counts in memory.
-  Disable counter so you can write totals and notes.

*design • function • innovation*



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

