



RPF05043 Solid Door Pharmacy & Vaccine Refrigerator*

Robustly constructed from Stainless Steel this refrigerator from the **Labcold** Advanced range is designed to provide close temperature control for temperature sensitive products.



The fan circulated interior has no internal 'step', making it easy to clean and stock and to allow for maximum air flow for even temperature distribution. The two shelves are fully adjustable, so it is possible to maximise the space available.

Insulated with Envirofoam, which is derived from plant rather than petrochemical oils, the refrigerant is also ammonia, HCFC and CFC free.

Fitted with a sturdy door lock, this Advanced refrigerator is ideal for busy environments such as busy laboratories or pharmacies.

Part Number	RPF05043
Capacity (Litres)	150
Style	Under Counter
Exterior Dimensions (HxWxD)	820 x 600 x 620mm
Interior Dimensions (HxWxD)	680 x 540 x 350mm
Energy Consumption	3.75 kWh/24h
Operating Temperature	+4°C
Temperature Set Point Range	+2°C to +8°C
Interior Construction	Stainless Steel
Exterior Construction	Stainless Steel
Shelves	2
Shelf dimensions (W x D)	530 x 325mm
Door Hinging	Right hand only
Door Style	Solid (Glass Door Model available RPF05043)
Internal Light	Glass Door Model only
Rear Castors	✓
Door Lock	✓
Access Port	Optional Extra
Digital Temperature Control	✓
Digital Temperature Display	✓
High Temperature alarm	✓
Low Temperature alarm	✓
Door Open alarm	✓
Auto defrost	✓
Fan Assisted Cooling	✓
Environmentally friendly insulation	✓
Eco Friendly Refrigerant	✓
Ammonia Free	✓
DH Green Book, RPS/ GPhC, NaTHNaC and RCVS Guidelines compliant	✓

* Glass door model available. Please order RPF05043.





Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

