



Basic Ward Range

Under Counter Fridge



FEATURES	KEY BENEFITS
Lock	Security is essential so our Lec Medical fridge comes complete with door lock and two keys, ensuring that only authorised individuals are able to access it.
Under counter size	Small enough to fit on a ward without being obtrusive.
12 month warranty	Our under counter fridge and freezer come complete with a 12 month warranty for added peace of mind.

SPECIFICATION TABLE

	Product Description		Basic Ward Fridge
 <p>Internal Drawer</p>	Product Dimensions (external)	mm	545 x 554 x 845mm (W x D x H)
	Capacity	litres	130
 <p>Lock</p>	Colour/Finish		White
	Reversible Door	Y/N	Yes
	Door Lock	Y/N	Yes
	Drawer Number		1 Plastic
	Shelf Number		3 Glass
	Handle Type		Recessed Handle
	Defrost Capability		Auto Defrost
	Annual Energy Consumption	kwh	117
	Thermostat		Adjustable
	Spark Free	Y/N	Yes
	HFC Free	Y/N	Yes
	CFC Free	Y/N	Yes





Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

