


NEW 1800mm wide Class 2 Biological Safety Cabinet with Visionaire™ Touch Screen Technology



Guardian MSCT1800 Specifications:

External Dimensions:	1800mm wide x 750mm deep x 2230mm high inc. base stand (1330mm without base stand) <small>(add 66mm height for outlet filter if fitted)</small>
Internal Dimensions:	1720mm wide x 630mm deep x 754mm high
Main Filter:	H14 HEPA filters – 99.997% efficient at 0.3µm
Primary Outlet Filter:	H14 HEPA filters – 99.997% efficient at 0.3µm
Secondary Outlet Filter:	OPTIONAL: H14 HEPA filter – 99.997% efficient at 0.3µm
Air Cleanliness:	ISO Class 4
Construction:	Epoxy coated Zintec Mild Steel with toughened glass glazing panels and removable 5-part 316 stainless steel work surface

**Monmouth
Scientific**





Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax: 01759 301143

sales@wolflabs.co.uk



Use the above details to contact us if this literature doesn't answer all your questions.

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

